

PIPEP™ BREATHING EXERCISER



- Latex – Free
- PVC Material
- Mask not included



Features and Benefits

- Respiratory exerciser system for inspiratory and expiratory breathing exercises with resistors for PEP or PIP.
- Equipped with one main connector, manometer connector and with 8 color coded resistors.
- Resistor holes' diameter : 1.5mm, 2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm, 4.5mm, 5.0mm
- Single patient use.
- Can be used with standard face mask



KOO MEDICAL EQUIPMENT

PIPEP™ BREATHING EXERCISER

Technical Characteristics

Product Number	Product Description	Product Size	Quantity Per Box	Description of contents
KM-803	PipeP	Standard	24	<ul style="list-style-type: none"> • Respiratory exercise system – in and out connector • 8 resistors • 1 manometer connector, OD 15m • In adapter size : OD 22mm • Out adapter size : OD 22mm ID 15mm



NOTE :

- Resistors can be used on inspiratory side or expiratory side to exercise different accessory muscle sets

Distributed by:



Latex free



KOO EUROPE SRL
 Via della Praia 12/D
 10090 - Buttigliera Alta (TORINO)
 ITALY
 Tel : +39 011 934 8192
 Fax : +39 011 934 8627
 Info@kooeuropa.com

KOO AMERICAS, INC
 301 East Wallace Kneeland Blvd
 Building 224 - 112
 Shelton, WA 98584
 Tel : +360 426 3669
 Fax : +360 358 1338
 andrea@koo-shanghai.com

KOO ASIA LTD.
 Unit 1, 16/F, Cable TV Tower,
 9 Hoi Shing Road, Tsuen Wan,
 N.T. HONG KONG
 Tel : +852 2110 4780
 Fax : +852 2110 4783
 sales@koo-asia.com